

**EZ-GO (P321) Group 1  
Super Light Mini RWD Power Wheelchair**

**Order Form** (Effective April 19, 2017)

For ease of ordering contact Inside Sales toll free at 1-800-963-7487, Fax: 239-574-2661



Order Taken by: \_\_\_\_\_

Company Name: \_\_\_\_\_ Account # \_\_\_\_\_

Phone: \_\_\_\_\_ Date: \_\_\_\_\_ P.O # \_\_\_\_\_

Ship to Name/Address: \_\_\_\_\_

**Special Note:**

DME suppliers are responsible for determining appropriate billing codes when submitting for insurance reimbursement. HCPCS codes should not be considered as legal advice.

**CHAIR TYPE PRICE HCPCS PRODUSPECIFICATIONS**

**EZ-GO Group 1**

P321-1S C Rear Wheel Drive Base( 18W X 16D ) \$3999.00 K0814

**FRONT RIGGINGS**

Footplate\* ..... No Charge

**JOYSTICK MOUNTING HARDWARE**

Swing Away Joystick Mount .....  
 Right Mount .....  Left Mount..... \$240.00 E1028

**BATTERY OPTIONS**

12V / 15A Battery (On Chair).....  
 \$55.00ea K0733

**ARM RESTS**

Flip Up, Adjustable-Height .....  
 \$150.00ea E0973

**SEATING OPTIONS (for P320-2S S)**

**OPTIONS & ACCESSORIES**

Oxygen Tank Holder..... \$138.00 E2208  
 Cane & Crutch Holder..... \$98.00 E2207  
 Walker Holder..... \$118.00 K0108  
 Foldable Rear Basket (Med.)..... \$72.00  
 Foldable Rear Basket (Lrg.)..... \$108.00  
 Cover..... \$180.00  
 Cup Holder\*\*..... \$24.00  
 Mounting Bracket..... \$32.00  
 Side Bag..... \$18.00

**SHROUD COLORS\***

Red..... Blue..... No Charge

<b>Dimensions</b>	
Length	34"
Width	21.5"
Height	33.5"
Ground Clearance	1.8"
<b>Seat Dimensions</b>	
Width	18"
Depth	16"
Back Height	15"
<b>Weight Capacity</b>	250-lbs
<b>Maximum Speed</b>	4 mph
<b>Range Up To</b>	9 miles
<b>Tires</b>	
Front	6" x 2"
Drive Wheel	9" x 3"
<b>Turning Radius</b>	23"
<b>Incline Capability</b>	6°
<b>Product Weight</b>	
Total Weight without Battery	79-lbs
Total Weight with Battery	99-lbs
Weight of Heaviest Piece	35-lbs
<b>Battery</b>	
Type	12V / 15A gel x 2
Weight	10-lbs x 2
<b>Warranty</b>	
Frame	5 year / Limited Lifetime
Electronics	1 year
Motor/Brakes	1 year
Batteries	6 montha

\*Standard with Unit

\*\*Mounting Bracket Required



Merits Health Products, Inc.,  
 730 NE 19<sup>th</sup> Place,  
 Cape Coral, FL 33909  
[www.meritshealth.com](http://www.meritshealth.com)

Quantity for above selection

**Features:**

- Light weight and easy to disassemble for handling, storage, and transportation
- Wireless electronic connector
- Ideal as personal secondary portable power wheelchair
- Compact size to easily fit into the trunk of most cars
- Solid Caster Wheels & Drive Wheels
- 1.5A off-board charger
- Lap Belt (Standard with chair)