

4-Wheel Stability with 3-Wheel Maneuverability

Jazzy[®]
ZERO TURN



Patented, intelligent-turning

iTurn
TECHNOLOGY™



Get a 4-wheel scooter with a 3-wheel turning radius with the **Jazzy® Zero Turn** and Pride's patented, intelligent-turning, iTurn Technology™. Navigate tight corners and small spaces effortlessly with a 38" turning radius. Illuminate your ride with bright LED lighting. Dual motors and CTS suspension ensure a smooth and comfortable drive.

PRIDE[®] MOBILITY SCOOTERS
LIVE YOUR BEST[®]

| pridemobility.com |





Specifications

Weight Capacity	300 lbs.
Maximum Speed ¹	Up to 4 mph
Ground Clearance ²	1.4" at motor, 1.5 at anti-tip, 3.5" (side deck)
Turning Radius ²	38.25"
Overall Length ²	45.75"
Overall Width ²	20"
Front Tires	2" x 7" solid
Rear Tires	2.5" x 9" solid
Range Per Charge ^{2,4}	Up to 10 miles
Total Weight Without Batteries	120 lbs.
Heaviest Piece When Disassembled	48 lbs. (rear section)
Standard Seating	Type: Foldable, compact, molded plastic Weight: 27 lbs. Material: Black vinyl Dimensions: Width: 17" Depth: 17"
Motor Output	24-volt DC motor, rear-wheel drive, two-motor
Brakes	Electronic, regenerative and electromechanical
Battery Requirements ^{3,5}	Size: (2) 22 AH Weight: Battery box 31.5 lbs. 13.5 lbs. each
Battery Charger	Off-board, 3.5 amp
Delta Tiller	Yes
Colors	

*Note color may appear more silver with hues of green in some lighting

Features

- The first 4-wheel indoor scooter
- Patented iTurn Technology™ provides 4-wheel stability with a 3-wheel turning radius
- Charging port built into the tiller for easy charging
- Contemporary new rim design makes the Jazzy® Zero Turn look as good as it drives
- CTS suspension for comfortable long range use
- Delta tiller, ergonomic wraparound handles enables easy-to-use left or right hand operation
- Dual motors for added traction in rough terrain
- Feather touch disassembly permits simple frame separation with only one hand
- Standard front basket
- Full LED lighting package for safe operation in low light settings
- Standard dual under seat storage bins
- USB charger built into the tiller for convenient charging of smart phones and portable devices

Accessories*

- Back mounted accessories
 - Rear basket
 - Oxygen holder
- Cell phone holder
- Cup holder
- Double rear mount, 3"
- Forearm crutch holder
- 50" lap belt
- Micro, black, weather cover
- Rear view mirror
- Saddle bag
- Safety flag
- Single rear mount, 3"
- Walker holder
- Wishbone crutch holder

*Only one rear accessory can be used at a time.

Option

- 22 AH battery box (additional spare)

Military Patches



- 1 Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition. These specifications can be subject to a variance of (+/- 10%).
- 2 Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or - 3%).
- 3 AGM or gel cell type required.
- 4 Tested in accordance with ANSI/RESNA, WC Vol2, section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.
- 5 Battery weight may vary based on manufacturer.

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Dimensions listed are for scooter only. Overall measurements will vary based on seating and accessory selections.

HCPCS Code: K0899

POWER MOBILITY DEVICE, NOT CODED BY DME PDAC
OR DOES NOT MEET CRITERIA

(US) 800-800-8586

(CAN) 888-570-1113

pridemobility.com



Pride[®]
Mobility Products Corp.