



**Titan X23  
Standard Powerchair**

HCPCS: K0823



INSET A



INSET B



INSET C



INSET D

- Height-adjustable swivel seat features fold-down backrest (Figure A)
- Includes interchangeable color panels in red and blue
- Seat back reclines and features adjustable headrest
- Adjustable-length controller mount (Figure B)
- Armrests are padded, adjustable and removable
- Available with 16", 18", 20" and 22" wide Captain's Seat
- Optional swing-away controller arm (Item #AA4100) sold separately
- Solid, flip-up footplate with adjustable height, depth and angle
- Positioning belt included
- Optional elevating leg rests (requires Item #AE2500 & Item #LK3JELR) sold separately
- Programmable controller
- Flat-free, non-marking tires (Figure C)
- Easy free wheel operation (Figure D)
- Anti-tip wheels

PRODUCTS	Description
TITANLB18CSX23	Titan X23 K0823 FWD 18 x 18 Capt.Seat, 1/cs
TITAN18CSX23	Titan X23 K0823 FWD 18 x 18 Capt.Seat, 1/cs
TITAN20CSX23	Titan X23 K0823 FWD 20 x 18 Capt.Seat, 1/cs
TITAN1616X23	Titan X23 K0823 FWD 16 x 16 Capt.Seat, 1/cs
TITAN1618X23	Titan X23 K0823 FWD 16 x 18 Capt.Seat, 1/cs
TITAN1816X23	Titan X23 K0823 FWD 18 x 16 Capt.Seat, 1/cs
TITAN2020X23	Titan X23 K0823 FWD 20 x 20 Capt.Seat, 1/cs
TITAN22CSX23	Titan X23 K0823 FWD 22 x 20 Capt.Seat, 1/cs

SPECIFICATIONS	Description
HCPCS	K0823
Weight Capacity	300 lbs.
Top Speed	4 mph
Maximum Range	15 miles
Turning Radius	25"
Climbing Angle	6°
Ground Clearance	2.5"
Dimensions	40" (L) x 24" (W)
Floor -to-Seat Height	20.5" - 23"
Controller	Dynamic Shark 40A
Motors	(2) 24V x 200W x 3900 rpm
Batteries (pair)	12V x 35AH (U1)

SPECIFICATIONS	Description
Charger	3A Offboard
Brakes	Electromagnetic
Caster Wheels	6" x 2" Solid
Drive Wheels	10" x 3" Flat-free
Base Weight	68 lbs.
Battery Weight (pair)	44 lbs.
Seat Weight	43 lbs.
Heaviest Piece Disassembled	55 lbs.
Total Weight	158 lbs.
Warranty on Frame	Limited Lifetime
Warranty on Electronic Controller/Drive Train Components	14 months
Warranty on Batteries	12 Months

ACCESSORIES	Description
AA4100	Swing-Away Joystick, 1/ea
AE2500	ELR Bracket, 1/ea
LK3JELR	ELR Set, 1/pr
AB1000	Arm Bag, 1/ea
AZ0060	Cup Holder, 1/ea