

# Vision Ultra

The Vision Ultra is a compact and versatile power wheelchair. With 6 wheels on the ground this unit has the stability and durability to withstand the test of time. This unit features 14" drive wheels, a 300 pound capacity and a seat-to-floor height starting at 16.5". The Vision Ultra can be ordered with 3 different seating systems to accommodate many users



# Vision Ultra

## Features

- Power tilt operated in 0°-55°
- Height-adjustable cantilever armrests
- Seat width and depth adjustable 16"-22" in 1" increments
- Large 14" drive wheels
- Tension-adjustable front and rear suspension

## Specification

Weight Capacity	300 lbs
Maximum Speed	5 mph
Per-Range Charge	12 miles (19 miles w/ 22NF)
Batteries	34Ah/U1
Ground Clearance	2.4"
Turning Radius	26"
Overall Length	36"
Overall Width	24"
Drive Wheels	14"
Casters	6" front, 6" rear
Suspension	Yes full
Seat Width	16"-22"
Seat Depth	16"-22"
Backrest Height	21"- 23" Tilt only, 20" Captain or pan
Seat-to-Floor Height	16.5"-19.5"
Warranty	18 months: Motor & Electronics 5 Year/Limited Lifetime Frame 1 Year Battery
Controller	P&G nVR2 and VR2-A2

## Options

- Captain Seat, Solid Seat Pan, or Tilt Seating
- Weather Cover
- Battery Upgrade: 22NF
- Swing- Away Joystick
- Articulating Power Footplate
- Trunk Supports (LxWxH)
  - Small: 5"x3"x1"
  - Med: 5"x4"x1"
  - Large: 6"x4"x1"
- Abductor (LxWxH)
  - Small: 4"x3.5"x3"
  - Large: 6"x4"x3"



## HCPCS

- Group 2 standard (Captain or pan seat) K0822/K0823
- Group 2 standard, single power: K0835

