

P O W E R C H A I R



Seat Colors
Grey



Body Colors
Red



Roll

Alanté
BY GOLDEN

The Alanté name continues to mean quality and value! The first thing that will catch your eye with the all-new Alanté is the unique shroud design that is stylish and functional. Featuring a full size captain's seat, the Alanté is a front wheel drive power base with our black, non-marking drive tires on high strength steel rims with bright chrome wheel covers.

Rise. Roll. Ride.

GOLDEN
technologies

YOUR LIFE IN MOTION

Alanté GP-207 Specifications

Model Number	GP207CC	GP207SS
Medicare Code	K0816	K0815
Weight Capacity	300 lbs	300 lbs
Drive Wheels	Front	Front
Maximum Speed	3 mph	3 mph
Operating Range ¹	17 miles	17 miles
Ground Clearance	1.75"	1.75"
Turning Radius	25"	25"
Type Batteries	2-U1	2-U1
Colors	Red	Red
Freewheel Mode	Yes	Yes
Electronic Speed Control	Yes	Yes
Electro-Mechanical Brakes	Yes	Yes
Charger	Off Board	Off Board
Controller	PG VR-2 50 Amp Remote	PG VR-2 50 Amp Remote
Dimensions:		
Length	39.5"	39.5"
Width	23.75"	23.75"
Height (ground to top of back of seat)	37" - 39"	41" - 43"
Height (ground to top of back of headrest)	44.5" - 48"	48.5" - 52"
Ground to Top of Seat	19" - 22"	17" - 20"
Number of Seat Height Adjustments	4	4
Size of Increments	1"	1"
Weight of Unit: (assembled)	163 lbs	160 lbs
Base	67 lbs	67 lbs
Seat w/Arms	46 lbs	43 lbs
Batteries	22 lbs ea.	22 lbs ea.
Tires:	Flat Free	Flat Free
Front	9"	9"
Casters	6"	6"
Front Anti-Tip Wheels	4"	4"
Standard Seat:	Captain's Seat	Captain's Van/Pan
Slide Seat	NA	NA
Back Height - no headrest	21.5"	23.5"
Back Height with headrest	26.5" - 30"	30.5" - 34"
Width x Depth	17.7" x 19.1"	17.7" x 19.1"
Color	Grey Vinyl	Grey Vinyl
Footrest	Yes	Yes
Height adjustable	Yes	Yes
Angle adjustable	Yes	Yes
Warranty:		
Frame	5 Years	5 Years
Drive Train	13 Month	13 Month
Electronics ²	13 Month	13 Month

¹. Battery range will vary due to rider weight, drive surface and terrain.

². Electronics warranty excludes batteries.

Features

- Low center of gravity
- Most stable unit in its class
- Excellent maneuverability
- 25" turning radius
- 300 lb weight capacity
- Upgraded user friendly motor connections
- Auto shut-off charger system
- Electromagnetic braking system
- Padded flip up arm standard
- 9" driving wheels
- Enlarged 4" anti-tipper wheels

Power Chair Accessories and Options

- Walker Holder
- Cane Holder
- Oxygen Tank Holder



Disassembles easily into 5 pieces!