

VivaLIFT![®]
POWER RECLINERS
By Pride[®] Mobility Products Corp.

THE MOST
COMFORTABLE
RECLINER
IN THE WORLD!

THE 4-MOTOR RECLINER

WITH
COUNTLESS
POSITION
COMBINATIONS







Shown in Saville
Wine Fabric

LEGACY




Embrace the serenity of a quiet corner in your home with a cozy, luxurious recliner. Take time for yourself with the Legacy Collection from **VivaLift!**[®]. Put your feet up and fully stretch out on the footrest extension. Enjoy a good novel or turn on the television. No matter what activity you choose, infinite, lay-flat positioning ensures you'll always find your most comfortable and supportive position. Plus, you'll never misplace your book or the remote with the convenient side storage pocket.

Pride Mobility[®] (US) 800-800-8586 | (CAN) 888-570-1113
pridemobility.com | prideliftchairs.com |    





LEGACY

COLLECTION	LEGACY v.2	LEGACY v.3
POSITION	Infinite Lay Flat Relaxed Reclined Upward Lift	Infinite Lay Flat Relaxed Reclined Upward Lift
WEIGHT CAPACITY	375 lbs.	375 lbs.
SEAT-TO-FLOOR	19.5"	20.5"
SEAT DEPTH	21.5"	22"
SEAT WIDTH	22"	23"
TOP OF BACK TO SEAT	25"	26"
FULLY PADDED CHAISE	Yes	Yes
FOOTREST EXTENSION	Standard	Standard
SEATING TYPE	Total Comfort	Total Comfort
CHAIR HEIGHT	42"	43.5"
DISTANCE FROM WALL	16"	16"
OVERALL WIDTH	35"	36.5"
CHAIR WEIGHT	153 lbs.	160 lbs.
WARRANTY	7-Yr. Pro-rated	7-Yr. Pro-rated

FUNCTIONS AND FEATURES:

-  POWER HEADREST
-  POWER LUMBAR
-  USB REMOTE
-  LITHIUM BATTERY BACKUP
-  FOOTREST EXTENSION
-  POWER BACKREST

POSITION GUIDE:

-  RELAXED
-  RECLINED
-  UPWARD LIFT
-  INFINITE LAY FLAT

STANDARD FABRIC: SAVILLE



WINE

GREY

BROWN

Note: Measurements may vary +/- 1". Specifications subject to change without prior notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.

 = Indicates shortest available lead time

Recommended HCPCS: E0627-Seat lift mechanism, electric, any type DME providers are responsible for determining the appropriate billing codes when submitting for insurance reimbursement. HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. Payer coding, coverage and bundling guidelines may apply.



Pride Mobility® (US) 800-800-8586 | (CAN) 888-570-1113

pridemobility.com | prideliftchairs.com |    