

**VivaLIFT!**<sup>®</sup>  
POWER RECLINERS  
By Pride<sup>®</sup> Mobility Products Corp.

THE MOST  
COMFORTABLE  
RECLINER  
IN THE WORLD!

## THE 4-MOTOR RECLINER

WITH  
COUNTLESS  
POSITION  
COMBINATIONS







Shown in Astro  
Brown Fabric

## TRANQUIL

There is no greater joy than reading a good book in a comfortable, cozy recliner. The Tranquil Collection by **VivaLift!**<sup>®</sup> Power Recliners combines quiet comfort with gorgeous fabrics, for a rich and upscale look. Get the support you need for your neck, shoulders and lower back as you dive into your favorite story. The Tranquil Collection exudes style and class with contrast stitching. It's the perfect addition to your home.









**Pride Mobility**<sup>®</sup> (US) 800-800-8586 | (CAN) 888-570-1113  
pridemobility.com | prideliftchairs.com |    

# TRANQUIL



COLLECTION	TRANQUIL v.2
POSITION	Infinite Lay Flat Relaxed Reclined Upward Lift
WEIGHT CAPACITY	375 lbs.
SEAT-TO-FLOOR	19.5"
SEAT DEPTH	20.5"
SEAT WIDTH	22"
TOP OF BACK TO SEAT	29"
FULLY PADDED CHAISE	Yes
FOOTREST EXTENSION	Standard
SEATING TYPE	Total Comfort
CHAIR HEIGHT	42"
DISTANCE FROM WALL	16"
OVERALL WIDTH	37.5"
CHAIR WEIGHT	145 lbs.
WARRANTY	7-Yr. Pro-rated




## FUNCTIONS AND FEATURES:

-  POWER HEADREST
-  POWER LUMBAR
-  USB REMOTE
-  LITHIUM BATTERY BACKUP
-  FOOTREST EXTENSION
-  POWER BACKREST

## POSITION GUIDE:

-  RELAXED
-  RECLINED
-  UPWARD LIFT
-  INFINITE LAY FLAT

## STANDARD FABRIC: ASTRO



GREY

BROWN

MUSHROOM

Note: Measurements may vary +/- 1". Specifications subject to change without prior notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.

 = Indicates shortest available lead time

Recommended HCPCS: E0627-Seat lift mechanism, electric, any type DME providers are responsible for determining the appropriate billing codes when submitting for insurance reimbursement. HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. Payer coding, coverage and bundling guidelines may apply.

**Pride Mobility**® (US) 800-800-8586 | (CAN) 888-570-1113

pridemobility.com | prideliftchairs.com |    