

# Pride® Lifts

*Mobility Made Easy.*

## Outlander DE

The Outlander DE is a dual-entry, three-rail scooter lift compatible with both 3- and 4-wheel small to mid-size scooters. The lift features an integrated automatic scooter lockdown arm designed to provide a secure self-locking mechanism for compatible mobility scooters. The simple hitched-based installation requires little to no modification to your vehicle or mobility device and allows for easy removal. The drive-on/drive-off platform enables loading from either end of the lift and takes the hassle out of loading your mobility scooter in a busy parking area.

### FEATURES AND BENEFITS

- 350 lbs. lifting capacity
- Compatible with Class II, III, or IV hitches
- LED puddle light for ease of use at night
- Toggle switch for ease of up and down use
- Integrated license plate bracket
- Phosphate undercoating for weather resistance
- LED indicator for On/Off
- Standard hardware kit
- Manual override with tool kit

BASE WEIGHT <sup>4</sup>	84 lbs.
MAXIMUM SCOOTER/PC OVERALL WEIGHT	350 lbs. <sup>2</sup>
MAXIMUM PRODUCT WIDTH	27"
OVERALL PLATFORM LENGTH	51" (42" usable surface)
MAXIMUM TRAVEL DISTANCE FROM THE GROUND	21"
MOTOR	12 Volt DC
WARRANTY	3-Year Transferable
HITCH CLASS <sup>3</sup>	II, III, IV
OPTIONS	<ul style="list-style-type: none"><li>• Battery pack w/ charger</li><li>• Swing-away adaptor</li><li>• Strap kit</li><li>• Scooter hold-down arm</li></ul>
VEHICLE PLATFORM	Full-size van, Minivan, Truck, Sedan, SUV and Crossover
LIFTS CURRENT PRODUCTS <sup>1,2</sup>	Victory® 9 3 & 4, Victory® 10 3 & 4, Victory® 10 LX w/CTS, Victory® Sport & Maxima 3 & 4 And equivalent competitor products <b>Except: Pursuit®, Pursuit® XL, Raptor</b>

1. Vehicle dependent.

2. Lifting weight does not include accessories, upgraded seating and/or hardware.

3. Hitch weight carrying capacity varies by vehicle. Call for details.

4. Does not include any accessories or hardware.

To find a compatible lift for your vehicle and mobility product, use our new

**LIFT FINDER**

[pridemobility.com/lift-finder](http://pridemobility.com/lift-finder)



MADE IN  
**USA**

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.