


Commander 450

Interior & Exterior Hoist Lift – 3-Axis



USER-FRIENDLY OPTION IDEALLY SUITED FOR PICKUP TRUCKS AND VANS

The hoist-style lift comes standard with a 450 lb lifting capacity, and boasts quick load and unload speeds. Rotation, coupled with some of the best lifting capacity in the industry make the Commander a great choice for heavier power wheelchairs and scooters.

This lift delivers maximum durability, with the ability to lift and lower your mobility device to the ground curbside at the rear of the truck bed or from the van's cargo area.

70"

Max Travel Distance from Ground

450

Pound Capacity

3

Powered Functions

COMMANDER 450 3-AXIS PRODUCT SPECIFICATIONS

Weight Capacity	450 lb (204 kg)
Vehicle Platform(s)	Full-Size Van, Truck
Accommodates	Heavy Power Chair, Scooter
Securement Option(s)	Docking Device
Power	12V DC
Max Travel Distance from Ground	70" (178 cm)
Warranty	3-Year Transferable



harmar.com/paragon | 800.833.0478