

Poly-Fly High Strength, Lightweight Wheelchair/ Flyweight Transport Chair Combo



FIGURE A



FIGURE B



FIGURE C



TR16

- Can be used as standard self-propelled wheelchair or a transport chair all in one
- Quick-release 24" wheels can be removed to transition self-propelled chair to a transport chair (Figure A - B)
- All-aluminum frame
- Comes with two sets of wheel locks: One for use with wheelchair, one for use with transport chair
- Back folds down with deluxe back release (Figure C)
- Nylon upholstery
- Comes standard with swing-away footrests with composite foot plates and heel loops
- 8" casters
- Comes with anti-tippers
- STDELR-TF - Swing-away Elevated Leg rest, Tool-free, (Optional front riggings)
- Weight Capacity: 250 lbs

PRODUCTS

PRODUCTS	Description
TR16	Silver Frame and Black Upholstery, 1/cs
TR18	Silver Frame and Black Upholstery, 1/cs
TR20	Silver Frame and Black Upholstery, 1/cs

SPECIFICATIONS

SPECIFICATIONS	Description
Overall Width with 24" Wheels (TR16)	22.375" (open) 24.5" (open with mag-wheels) 12" (closed)
Overall Width with 24" Wheels (TR18)	23.625" (open) 26.5" (open with mag-wheels) 12" (closed)
Overall Width with 24" Wheels (TR20)	26.25" (open) 29" (open with mag-wheels) 12" (closed)
Overall Height w/ 24" Wheels	37"
Seat-To-Floor Height	19.5"
Back Height	16"
Seat Depth	16"
Overall Length w/ Rigging	42"
Weight w/ 24" Wheels (TR16 TR18 TR20)	33 lbs. 34 lbs 36 lbs.
Carton Shipping Weight w/ 24" Wheels (TR16 TR18 TR20)	38 lbs. 39 lbs. 41 lbs.
Weight Limit	250 lbs.
Warranty	Limited Lifetime