



Aluminum Transport Chair



ATC17-BK



ATC17-SL



ATC17-RD

- Weighs only 19 lbs
- Aluminum frame is lightweight and strong
- Composite 8" wheels are lightweight and maintenance free
- Both 17" and 19" come with swing-away footrests
- Comes with seat belt for added safety
- Padded armrests provide additional comfort
- Nylon upholstery is comfortable, lightweight and easy-to-clean Back folds down for storage and transport
- Aluminum casted back-release hinge allows back to fold down for easy storage and transport

PRODUCTS

Description

ATC17-BK	Black Frame and Black Upholstery, Transport Chair17 Black, 1/cs
ATC17-BL	Blue Frame and Black Upholstery, Transport Chair17 Blue, 1/cs
ATC17-RD	Red Frame and Black Upholstery, Transport Chair Alum.17Red, 1/cs
ATC17-SL	Silver Frame and Black Upholstery, Transport Chair17Silver, 1/cs
ATC19-BK	Black Frame and Black Upholstery, Transport Chair Blk 19 Black, 1/cs
ATC19-BL	Blue Frame and Black Upholstery, Transport Chair Alum.19in Blue, 1/cs
ATC19-RD	Red Frame and Black Upholstery, Transport Chair Alum.19Red, 1/cs
ATC19-SL	Silver Frame and Black Upholstery, Transport Chair Alum.19 Silver, 1/cs

SPECIFICATIONS

Description

Overall Width (ATC17)	20" (open) 9" (closed)
Overall Width (ATC19)	22" (open) 9" (closed)
Seat-To-Floor Height	19"
Seat Dimensions (ATC17)	17" (W) x 16" (D)
Seat Dimensions (ATC19)	19" (W) x 16" (D)
Weight	19 lbs.
Overall Length w/ Riggings	33"
Overall Height	39.5"
Carton Shipping Weight	24 lbs.
Weight Limit	300 lbs.
Warranty	Limited Lifetime

ACCESSORIES

Description

STDS1040S	Universal Cup Holder, 1/ea
-----------	----------------------------