

Deluxe Go-Kart Steel Transport Chair (Chrome)




FIGURE A

FIGURE B

- Weighs 26 lbs
- Carbon steel frame with double-coated chrome provides an attractive, chip-proof, easy-to-maintain finish (Figure A)
- Composite 8" wheels are lightweight and maintenance free
- Both 17" and 19" come with swing-away footrests
- Padded armrests provide additional comfort
- Nylon upholstery is comfortable, lightweight and easy-to-clean
- Back folds down for storage and transport
- Aluminum casted back-release hinge allows back to fold down for easy storage and transport
- Comes with seat belt for added safety
- Comes with carry pocket on backrest
- Comes with upholstery handles to fold chair
- Standard rear wheel locks
- Convenient cup holder can easily attach to chair
- Sold separately (Model #STDS1040S)

PRODUCTS

Description

TR19	Chrome Frame, Black Upholstery, 1/cs
------	--------------------------------------

SPECIFICATIONS

Description

Overall Width (TR17)	20" (open) 9" (closed)
Overall Width (TR19)	22" (open) 9" (closed)
Seat-To-Floor Height	19"
Back Height	18"
Seat Dimensions (TR17)	17" (W) x 16" (D)
Seat Dimensions (TR19)	19" (W) x 16" (D)
Overall Length with Riggings	33"
Weight	26 lbs.
Carton Shipping Weight	31.5 lbs.
Weight Limit	300 lbs.
Warranty	Limited Lifetime

ACCESSORIES

Description

STDS1040S	Universal Cup Holder, 1/ea
-----------	----------------------------

PARTS

Description

STDS2351BLK	back upholstery for s/n with A, 1/bx
TRHG	Handgrip, 1/bx
TRB-1	Back Release Lever, 1/ea
TRH	folding hinge, 1/ea

PARTS	Description
STDS002001	arm pad bolts, 1/ea
STDS2350BLK	seat upholstery, 1/bx
STDS4011	rubber tip, 1/ea
STDS00200R	right side brake, 1/ea
STDS00200L	left side brake, 1/ea
TR-AXLE	8" Rear Wheel Axle, 1/ea
STDS002050	8" Rear Wheel Bearing, 1/ea
TRW-1	rear wheel, 1/ea
STDS1002-1	Caster axle bolt, 1/ea
STDS1011	Caster Bearing, 1/ea
STDS1066	Dust Cap, 1/ea
STDS4020	tipping sleeve, 1/ea
STDS1002	Standard 8" Caster, 1/ea
STDS1060	Fork Stem Bearing, 1/ea
STDS001045	Fork Stem Nut, 1/ea
STDS001000	Axle Cap, 1/ea
FW-1001	front seat rail guide, 1/ea
STDS1070	Seat Rail Guide, 1/ea
STDELR-2L	release lever-left, 1/ea
STDELR-2R	release lever- Right, 1/ea
STDS00202R	Footrest hanger (right), 1/ea
STDS00202L	Footrest hanger (left), 1/ea
STDS001051	Footrest adjustment bolt, 1/ea
STDS1026	Footrest Bumper, 1/ea
FWFP-R	Footplate- right, 1/ea
STDS001061	cross brace bolt, 1/ea
STDS002044	arm rest side panel, 1/ea
STDS5B	arm release button, 1/ea
STDS1M2017	for s/n with M, 1/ea
TR-17B	back upholstery for s/n with A, 1/ea
STDS5SKNB	Desk arm pad, 1/ea
STDS5FLNB	Arm Pad- Full, 1/ea
TR-17S	seat upholstery for s/n with A, 1/ea
STDS850-PB	Seat Belt, 1/ea
FWFP-L	Footplate- left, 1/ea
STDS002017	Cross Frame/TC, 1/ea
STDS1072-8	Seat Upholstery Scews (8), 1/ea
STDS5D5504	Back Upholstery Screws (8), 1/pkg