



Blue Streak Wheelchair Single Axle

HCPCS: K0001

- Blue Powder-coated frame, black cross brace
- Flip-back, padded desk-length arm
- Nylon upholstery is durable, lightweight, attractive and easy to clean
- Solid rubber tires mounted on composite wheels provide durability and low maintenance
- Comes with push-to-lock wheel locks
- Comes with calf strap
- Comes with swing-away footrests or elevating leg rests
- Weight Capacity: 250 lbs

PRODUCTS

Description

BLS18FBD-SF	BlueStreak18D/L FlipBackArms SwgFt, 1/cs
BLS18FBD-ELR	BlueStreak18D/L FlipBackArms Elev, 1/cs
BLS16FBD-SF	BlueStreak16D/L FlipBackArms SwgFt, 1/cs
BLS16FBD-ELR	BlueStreak16D/L FlipBackArms Elev, 1/cs
BLS20FBD-SF	BlueStreak 20D/FlipBackArms SwgFt, 1/cs
BLS20FBD-ELR	BlueStreak 20D/FlipBackArms ElevLeg, 1/cs

SPECIFICATIONS

Description

Width Open with 18" Seat	24"
Width Open with 20" Seat	26"
Width Closed	12.5"
Height	36"
Seat-To-floor Height	19.5"
Back Height	16"
Seat Dimensions with 18" Seat	18" (W) x 16" (D)
Seat Dimensions with 20" Seat	20" (W) x 16" (D)
Length with Riggings	42"
Product Weight	41.25 lbs.
Carton Shipping Weight	47 lbs.
Weight Limit	250 lbs.
Warranty	Limited 3 Year

ACCESSORIES

Description

STDS806	Anti-Fold Bar, 1/ea
STDS833N	Anti-Tippers With Wheels Black, 1/pr
STDS831	Heel Loops, 1/bx
WASR	Limb Support, Right, 1/ea
WASL	Limb Support, Left, 1/ea
STDS850	Seat Belts, Auto-Clasp Type Closure, 1/ea

ACCESSORIES	Description
STDS851	Seat Belts, Hook and Loop Fastener Type Closure, 1/ea
STDS855	Seat Belts, Auto-Clasp Type Closure, Bariatric, 1/ea
STDS856	Seat Belts, Hook and Loop Fastener Type Closure, Bariatric, 1/ea
STDS801GT	Wheel Lock Extension, 1/ea
STDS1034	Wheelchair Cane/Crutch Holder, 1/ea
PARTS	Description
STDS806	Anti-Fold Bar, 1/ea
STDS833	Anti-Tippers With Wheels, 1/pr
STDS833N	Anti-Tippers With Wheels Black, 1/pr
STDS831	Heel Loops, 1/bx
WASR	Limb Support, Right, 1/ea
WASL	Limb Support, Left, 1/ea
STDS850	Seat Belts, Auto-Clasp Type Closure, 1/ea
STDS851	Seat Belts, Hook and Loop Fastener Type Closure, 1/ea
STDS855	Seat Belts, Auto-Clasp Type Closure, Bariatric, 1/ea
STDS856	Seat Belts, Hook and Loop Fastener Type Closure, Bariatric, 1/ea
STDS801GT	Wheel Lock Extension, 1/ea
STDS1034	Wheelchair Cane/Crutch Holder, 1/ea
STDS4S2209	1/ea
STDS4S2208	1/ea
STDS4S2210	1/ea
TRHG	1/bx
STDS3J2424	1/ea
STDS8M481R	1/ea
STDS3J2400	1/ea
STDS 3J24131	
STDS 3J2421	
STDS 4S2225	
STDS 4S221R	
STDS 4S221L	
STDS8M4603	1/ea
STDS5N2202	2/pkg
STDS3J2416	1/ea
STDS3J2418	1/ea
STDS3J2420	1/ea
STDS4S22X8	1/ea
STDS4S2206	1/ea
STDS4S2207	1/ea
STDS5N2231	2/pkg
STDS3J 2432	
STDS3J4533	1/ea
STDS8M4601	1/ea
STDS4S22SF	4/pkg
STDS5N2201	1/pkg
STDS5N2210	1/ea
STDS1070	1/ea
STDSDSKVB	1/ea
STDS5N2219	1/pkg
STDS5N22DL	1/ea
STDS5N22DR	1/ea
STDS4S2211	1/ea
STDS4S2212	1/ea
LK3JELR	1/pr

- Comes with swing-away footrests or elevating leg rests

STDS004609-9	1/bag
STDS1026	1/ea
STDS1072-8	1/ea
STDS3J2413	1/ea
STDS3J2421	1/ea
STDS3J2432	1/ea
STDS4S0906-1	
STDS4S1069	1/ea
STDS4S221L	
STDS4S221R	
STDS4S2225	
STDS4S222X	
STDS4S22X6	
STDS5N2235	1/ea
STDS8M4925	1/ea
WB 8036	1/ea