



## Silver Sport 1 Wheelchair Single Axle

HCPCS: K0001




FIGURE A



INSET B



INSET C



INSET D



INSET E



INSET F

- Powder coated silver vein steel frame (Figure A)
- Padded armrests
- Nylon upholstery is durable, lightweight, attractive and easy-to-clean
- Urethane tires mounted on composite wheels provide durability and low maintenance
- Comes with push-to-lock wheel locks
- Carry pocket on backrest provides additional convenience
- Comes with swing-away footrests
- Folding back-canes for transport and storage
- Weight Capacity: 250 lbs

### PRODUCTS

#### Description

SSP118FA-SF	Silver Vein, 1/cs
-------------	-------------------

### SPECIFICATIONS

#### Description

Overall Width (open)	24"
Overall Width (closed)	11"
Overall Height	36"
Seat-To-Floor Height	19.5"
Back Height	16"
Seat Dimensions	18" (W) x 15.75" (D)
Overall Length (with riggings)	43.5"
Weight	37.6 lbs.
Carton Shipping Weight	42.8 lbs.
Weight Limit	250 lbs.
Warranty	Limited 3 Year
Overall Height (folded)	30"

### ACCESSORIES

#### Description

STDS806	Anti-Fold Bar, 1/ea
STDS833N	Anti-Tippers With Wheels Black, 1/pr
STDS818	Anti-Tippers without Wheels, Chrome, 1/pr
STDS829	Anti-Tippers With Wheels, Chrome, 1/pr
STDS831	Heel Loops, 1/bx
WASR	Limb Support, Right, 1/ea
WASL	Limb Support, Left, 1/ea
STDS834	Overhead Anti-Theft Device Single Pole, 1/ea
STDS850	Seat Belts, Auto-Clasp Type Closure, 1/ea
STDS851	Seat Belts, Hook and Loop Fastener Type Closure, 1/ea
STDS855	Seat Belts, Auto-Clasp Type Closure, Bariatric, 1/ea

## ACCESSORIES

## Description

STDS856	Seat Belts, Hook and Loop Fastener Type Closure, Bariatric, 1/ea
STDS820	Telescoping IV Pole Attachment, 1/ea
STDS862	Universal Forward Stabilizer, 1/pr
STDS860	Universal Forward Stabilizer, 1/pr
STDS6008-1	Oxygen Cylinder Carry Bag for Wheelchair, 1/ea
STDS801GT	Wheel Lock Extension, 1/ea
STDS1034	Wheelchair Cane/Crutch Holder, 1/ea
LK3JELR	ELR Set, 1/pr
STDELR-AL	Articulating Legrest For Manual Wheelchairs, 1/cs
STDSNUT	Tinnerman Nuts, 1/pr
STDS5050	Easy-To-Clean Surface. Accommodates Wheelchairs, 1/ea
STDS1040S	Universal Cup Holder, 1/ea
STDS6005-1	Carry Bag,W/C,Back,15x15,1/ea, 1/cs
14300	Tension Adjustable Back Cushion, 1/ea

## PARTS

## Description

STDS806	Anti-Fold Bar, 1/ea
STDS833N	Anti-Tippers With Wheels Black, 1/pr
STDS818	Anti-Tippers without Wheels, Chrome, 1/pr
STDS829	Anti-Tippers With Wheels, Chrome, 1/pr
STDS835	Chart Carry Pocket, 1/ea
STDS831	Heel Loops, 1/bx
WASR	Limb Support, Right, 1/ea
WASL	Limb Support, Left, 1/ea
STDS823	Overhead Anti-Theft Device, 1/ea
STDS834	Overhead Anti-Theft Device Single Pole, 1/ea
STDS850	Seat Belts, Auto-Clasp Type Closure, 1/ea
STDS851	Seat Belts, Hook and Loop Fastener Type Closure, 1/ea
STDS855	Seat Belts, Auto-Clasp Type Closure, Bariatric, 1/ea
STDS856	Seat Belts, Hook and Loop Fastener Type Closure, Bariatric, 1/ea
STDS820	Telescoping IV Pole Attachment, 1/ea
STDS860	Universal Forward Stabilizer, 1/pr
STDS803	Universal Oxygen "E" Cylinder Holder, 1/ea
STDS804	Universal Oxygen "E" Cylinder Holder with IV Pole Attachment, 1/ea
STDS6008-1	Oxygen Cylinder Carry Bag for Wheelchair, 1/ea
STDS801GT	Wheel Lock Extension, 1/ea
STDS1034	Wheelchair Cane/Crutch Holder, 1/ea
LK3JELR	ELR Set, 1/pr
STDELR-AL	Articulating Legrest For Manual Wheelchairs, 1/cs
STDSNUT	Tinnerman Nuts, 1/pr
STDS5050	Easy-To-Clean Surface. Accommodates Wheelchairs, 1/ea
STDS1040S	Universal Cup Holder, 1/ea
STDS6005-1	1/cs
14300	Tension Adjustable Back Cushion, 1/ea