



Viper Plus GT Wheelchair
Deluxe High-Strength, Lightweight, Dual-Axle

HCPCS: K0004



FIGURE A



FIGURE B



FIGURE C

- Built-in seat rail extensions and extendable upholstery easily adjusts seat depth from 16" to 18" (Patent pending) (Figure A - C)
- All-Aluminum frame
- Dual Axle allows for multiple seat-to-floor height positions
- Adjustable height, flip-back, removable, padded arms
- Weighs under 35 lbs in basic configuration (excluding front rigging)
- Over-center cross brace eliminates need for seat guides Works with seating and positioning accessories
- Adjustable angle back easily adjusts from 5° – 20°
- Adjustable height back adjusts from 17" – 19" in 1" increments, adaptable for individual patient needs
- Durable, heavy-gauge flame-retardant nylon upholstery with hook and loop fastener flaps
- Quick release, 24" composite Mag-style wheels are lightweight and maintenance free
- Precision sealed wheel bearings in front and rear ensure long-lasting performance and reliability
- Adjustable angle caster forks have three height adjustments and angle adjustability
- Tool-free adjustable swing-away footrests also feature easy-to-access release handles
- Urethane rear tires and casters offer superior performance, smoother ride and are lightweight
- Adjustable anti-tipper with flip-up wheels Never needs to be removed
- Comes with push-to-lock wheel locks
- Patent Pending

PRODUCTS

Description

PLA416FBFAARAD-SF	ViperPlsGT 16FlpBk.AdjHtFullArm SF, 1/cs
PLA416FBFAARAD-ELR	ViperPlsGT 16FlpBkAdjHtFullArm ELR, 1/cs
PLA416FBDAAARAD-SF	ViperPlsGT 16FlpBk.AdjHt.DskArm SF, 1/cs
PLA416FBDAAARAD-ELR	ViperPlsGT 16FlpBk.AdjHtDskArm ELR, 1/cs
PLA418FBDAAARAD-SF	ViperPlsGT 18FlpBk.AdjHt.DskArm SF, 1/cs
PLA418FBDAAARAD-ELR	ViperPlsGT 18FlpBkAdjHtDskArm EL, 1/cs
PLA418FBFAARAD-SF	ViperPlsGT 18FlpBk.AdjHtFullArm SF, 1/cs
PLA418FBFAARAD-ELR	ViperPlsGT 18FlpBkAdjHtFullArm ELR, 1/cs
PLA420FBDAAARAD-SF	ViperPlsGT 20FlpBkAdjHtDetDskArmSF, 1/cs
PLA420FBDAAARAD-ELR	ViperPlsGT20FlpBkAdjHtDetDskArmELR, 1/cs
PLA420FBFAARAD-SF	ViperPlsGT20FlpBkAdjHtDetFullArmSF, 1/cs
PLA420FBFAARAD-ELR	ViperPlsGT20FlpBkAdjHtDtFullArmELR, 1/cs
PLA422FBDAAAR-SF	ViperPlsGT 22FlpBkAdjHtDetDskArmSF, 1/cs
PLA422FBDAAAR-ELR	ViperPlsGT22FlpBkAdjHtDetDskArmELR, 1/cs
PLA422FBFAAR-SF	ViperPlsGT22FlpBkAdjHtDetFullArmSF, 1/cs
PLA422FBFAAR-ELR	ViperPlsGT22FlpBkAdjHtDtFullArmELR, 1/cs

SPECIFICATIONS

Description

Width Open with 16" Seat	24"
--------------------------	-----

