



Cougar Wheelchair
Ultralight Aluminum Wheelchair

HCPCS: K0005




FIGURE A



FIGURE B



FIGURE C

- Ultra-light, weighs less than 27 lbs (K0005)
- Forward folding back to reduce size
- Tool-free back height adjustable from 17"-20 ½"
- Seat depth adjustable from 16"-18"
- Multi-position back can adjust 15°
- Rigidity bar provides support
- Quick-release rear wheels (Figure A)
- Semi-Pneumatic tires on a 24" spoke wheel provide a smooth, true roll
- 3" of depth adjustment in the overall wheelbase
- Front casters are angle adjustable and height adjustable in three positions
- Arm rest can be flipped allowing 2" width adjustment
- Nylon upholstery is breathable, lightweight and comfortable
- Back folds forward for storage (Figure B - C)
- Adjustable aluminum wheel locks
- Swing-away footrests swing in and out
- Lightweight, composite foot plates with heel loops come standard
- Optional adjustable height arm (sold separately)

PRODUCTS

Description

AK516ADA-AELR	Cougar UltraLt W/C 16Det, 1/cs
AK516ADA-ASF	Cougar UltraLt W/C16Det TA, 1/cs
AK518ADA-AELR	Cougar UltraLt W/C 18Det, 1/cs
AK518ADA-ASF	Cougar UltraLt W/C 18Det, 1/cs

SPECIFICATIONS

Description

Width Open with 16" Seat	24"
Width Open with 18" Seat	26"
Width Closed	12"
Height	37"
Seat-To-Floor Height	17.5" - 19.5"
Seat Depth	16"
Length with Riggings	44"
Product Weight	27 lbs.
Carton Shipping Weight	32 lbs.
Weight Limit	250 lbs.
Warranty	Limited Lifetime

ACCESSORIES

Description

STDS601	Wheel Lock Extension, 1/ea
---------	----------------------------

ACCESSORIES

Description

AK5ADJARMKIT-L	Adjustable Height Arm, Left, 1/ea
AK5ADJARMKIT-R	Adjustable Height Arm, Right, 1/ea
STDS806	Anti-Fold Bar, 1/ea
STDS2A4326	With Wheels, 1/pr
STDS834	Overhead Anti-Theft Device Single Pole, 1/ea
STDS850	Seat Belts, Auto-Clasp Type Closure, 1/ea
STDS851	Seat Belts, Hook and Loop Fastener Type Closure, 1/ea
STDS820	Telescoping IV Pole Attachment, 1/ea
STDS860	Universal Forward Stabilizer, 1/pr
STDS1034	Wheelchair Cane/Crutch Holder, 1/ea
STDS831	Heel Loops, 1/bx

PARTS

Description

STDS601	Wheel Lock Extension, 1/ea
AK5ADJARMKIT-L	1/ea
AK5ADJARMKIT-R	1/ea
STDS806	Anti-Fold Bar, 1/ea
STDS2A4326	With Wheels, 1/pr
STDS835	Chart Carry Pocket, 1/ea
STDS831	Heel Loops, 1/bx
STDS823	Overhead Anti-Theft Device, 1/ea
STDS821	Overhead Anti-Theft Device With IV Hooks, 1/ea
STDS834	Overhead Anti-Theft Device Single Pole, 1/ea
STDS850	Seat Belts, Auto-Clasp Type Closure, 1/ea
STDS851	Seat Belts, Hook and Loop Fastener Type Closure, 1/ea
STDS820	Telescoping IV Pole Attachment, 1/ea
STDS860	Universal Forward Stabilizer, 1/pr
STDS1034	Wheelchair Cane/Crutch Holder, 1/ea
STDS2A4301	
STDS2A431R	
STDS2A431L	
STDS2A4306	
STDS2A4308	
STDS2A4335	
STDS2A43X8	
STDS2A432X	
STDS2A4309	
STDS2A4333	
STDS2A4304	
STDS2A4310	
STDS2A4313	
STDS2A4314	
STDS1M4047	
STDS2A4321	
STDS2A4322	
STDS2A432R	
STDS2A432L	
STDS4058	
STDS2A4302	
STDS2A4323	
STDS2A4324	
STDS2A4326	
STDS2A4330	
STDS2A4300	