

AL500 / AL500HD

Universal Power Chair Lift



DESIGNED FOR POWER CHAIRS

American-made AL500 and AL500HD Universal Power Chair Lifts are among the world's most popular outside platform lifts. Featuring Harmar's Stars-N-Stripes aluminum deck design, dedicated to our veterans and commemorating Harmar's heritage as an American company.

The AL500 and AL500HD are easy to operate and built to last for years of trouble-free use. The universal design safely transports the power chair without making chair modifications.

- Easy, one-switch operation
- License plate light
- Key switch and courtesy light
- 4 - Self-tensioning Q'Straint retractors, no chair modifications required
- Manual crank backup and license plate mount included
- Swing-Away option available

EASY OPERATION, BUILT TO LAST

AL500 / AL500HD



I am a caretaker for my father and this lift made it possible for me to take my dad to go out to eat and meet with friends and family. He now feels a sense of freedom and enjoys trips big and small.

75

Pound Install Weight

350

Pound Capacity

77

Pound Install Weight

400

Pound Capacity

AL500 Product Specifications

| | |
|------------------|---|
| Accommodates | Power Wheelchairs |
| Platform | 26" x 39" aluminum deck |
| Lifting Capacity | 350 lb |
| Installed Weight | 75 lb |
| Securement | Self-tensioning Q'Straint retractors |
| Hitch | II or greater |
| Options | Battery pack, Swing-Away, Power Wheelchair Covers |
| Warranty | 3-Year Transferable |

AL500HD Product Specifications

| | |
|------------------|---|
| Accommodates | Power Wheelchairs |
| Platform | 26" x 39" aluminum deck |
| Lifting Capacity | 400 lb |
| Installed Weight | 77 lb |
| Securement | Self-tensioning Q'Straint retractors |
| Hitch | III or greater |
| Options | Battery pack, Swing-Away, Power Wheelchair Covers |
| Warranty | 3-Year Transferable |

AL500 AL500HD / DA110918



Harmar.com | 800.833.0478

